






## HYPODERMIC SYRINGES IN 3 PARTS WITH AND WITHOUT NEEDLE

High-quality range of medical devices ensuring an optimal safety. The range of Hypodermic syringe INJ/LIGHT have been manufactured in compliance with the most important international standards in conformity with European Directive 93/42/EEC and with UNI EN ISO 9001/2000 - UNI EN ISO 13485/2004.


The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack. Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.



SYRINGE WITHOUT NEEDLE		
GIMA code	Capacity	Minimum order
	LUER LOCK CONE	
23808	10 ml	box of 100
23809	20 ml	
23810	30 ml	
23811	60 ml	
		
CENTRIC LUER CONE		
23813	3 ml	box of 100
23814	5 ml	
23815	10 ml	
ECCENTRIC LUER CONE		
23816	10 ml	box of 100
23817	20 ml	
23819	60 ml	
		
CATHETER CONE		
23820	60 ml	box of 20
23821	100 ml	box of 25
		

SYRINGE WITH NEEDLE			
GIMA code	Capacity	Assembled needle	Minimum order
	CENTRIC LUER CONE		
23822	1 ml - tubercolin	26G - 0.45x12.7 mm	box of 100
23824	2.5 ml*	21G - 0.8x38 mm	
23825	2.5 ml*	22G - 0.7x32 mm	box of 100
23827	5 ml	21G - 0.8x38 mm	
23828	5 ml	22G - 0.7x32 mm	box of 100
23830	10 ml	21G - 0.8x38 mm	
 <p>*Extended scale: 3 ml</p>			
ECCENTRIC LUER CONE			
23833	10 ml	21G - 0.8x38 mm	box of 100
23834	20 ml	21G - 0.8x38 mm	
			


Needle in AISI 304 stainless steel triple bevelling



## HYPODERMIC NEEDLES AND SCALP-VEIN SETS - STERILE

A line of hypodermic needles and scalp vein sets ensuring minimal trauma during penetration thanks to the triple bevelling. Each element of the needle is treated for painless injection.

The needles are in stainless steel AISI 304 with triple bevelling and totally non-toxic




HYPODERMIC NEEDLES				
Standard	Safety	Needle	Colour	Minimum order
23750	23780	18G - 1.2x38 mm	Pink	box of 100 (standard needles)
23751	23781	19G - 1.1x38 mm	Cream	
23753	-	20G - 0.9x38 mm	Yellow	
23755	23785	21G - 0.8x38 mm	Green	
23757	-	22G - 0.7x32 mm	Black	box of 50 (safety needles)
23758	23788	23G - 0.6x25 mm	Turquoise	
23760	-	25G - 0.5x16 mm	Orange	
23761	-	26G - 0.45x12.7 mm	Brown	
23762	-	27G - 0.4x12.7 mm	Grey	

and non-pyrogenic. Needles have Luer-slip or Luer-lock. The scalp vein set needles are provided with soft flexible grained wings. The tube is made of 30 cm PVC with luer-lock. Individually packaged in blister: GB, FR, IT, DE, ES, PT, GR, RU (Scalp-vein sets) GB, IT, ES, PT (Hypodermic needles).



SCALP - VEIN SETS				
Standard	Safety	Needle	Colour	Minimum order
23766	-	19G - 1.1x19 mm	Cream	box of 100
23767	-	20G - 0.9x19 mm	Yellow	
23768	23774	21G - 0.8x19 mm	Green	box of 100
23769	-	22G - 0.7x19 mm	Black	
23770	23776	23G - 0.6x19 mm	Turquoise	box of 100
23771	-	25G - 0.5x19 mm	Orange	

## 5 ml SYRINGE - RETAIL BOX

• 23795 SYRINGE 5 ml  
With needle 23G 0.60x25 mm.  
Sterile, packed in single blister.  
Box of 10 syringes.   
Minimum order:  
carton of 12 boxes of 10 syringes



## 5 ml SAFETY SYRINGE - RETAIL

• 23796 SAFETY SYRINGE 5 ml  
With needle 23G 0.60x25 mm.  
Sterile, packed in single blister.  
Box of 5 syringes.  
Minimum order:  
carton of 48 boxes of 5 syringes

**STOP ACCIDENTAL NEEDLE STICKS**

