

Projet 101 S2 Ear Irrigator



The Projet Ear Irrigator has been updated and improved. The new design increases the safety, comfort and effectiveness of ear canal cleaning by utilising its gentle, pulsating cyclone wash.

Every design aspect has been repeatedly tested clinically and readjusted until the best results have been achieved.

To enhance performance and cleaning efficiency, revised features include:

- Very simple to operate
- Smart protection with auto shut off
- Quiet operation promoting noise reduction •
- Improved pulsating frequency

Ear Irrigation kit consists of:

- Irrigator machine
- Water dispensing handle •
- Water reservoir
- 10 ear irrigation tips
- Mains transformer
- Foot switch

Specifications

- Power: 240V, 60Hz
- Static pressure: 0.25 - 3 bar maximum •
- Working time: 2 - 3 mins for 500cc
- Continuous usage: 5 mins
- Reservoir: 500ml
- Pressure: 860 - 1060 hpa

Code	Description
190.60.000/S2	Projet S2 Ear Irrigator
190.60.100	Projet S2 Ear Tips

Projet