

VITALOGRAPH myPEF PEAK FLOW METER

For Adults & Children



VITALOGRAPH myPEF


Peak Flow Meter

The effective management of asthma relies on early identification of the threat of potential exacerbations. The ability to respond promptly depends on close monitoring of peak flow to help protect quality of life and optimise the efficacy of medication.

Vitalograph® has been providing simple and effective peak flow meters for use in healthcare environments and the home for over 50 years. The Vitalograph myPEF™ Peak Flow Meter is easy to use and highly accurate. The attractive, ergonomic design and the simple colour coded scale help to support patient compliance and effective monitoring.



- Suitable for adults and children
- For use in clinics, general practice, sports medicine and at home
- Highly accurate
- Lightweight and portable
- Ergonomically designed with comfortable integral mouthpiece
- Features adjustable 'traffic light' colour coded scale and action plan
- Conforms to all relevant standards for peak flow meters
- 50-800L/min. measuring range
- Suitable for multi patient use, for training purposes, when used with SafeTway™ mouthpieces
- Includes a robust and hygienic protective case
- Branding options available



Authorised Distributor



Ordering Information

43900 myPEF Peak Flow Meter (1)

Actual product may vary from shown branding

Accessories*

20242 SafeTway Mouthpiece with one way valve (200)

29805 myPEF Adapter: for use with SafeTway Mouthpieces (96)

* Accessories only required if myPEF is used on multiple patients

Vitalograph Ltd, Intl., Maids Moreton, Buckingham, MK18 1SW, England T: +44 1280 827120 F: +44 1280 823302 E: sales@vitalograph.co.uk www.vitalograph.eu

Vitalograph®, SafeTway® and myPEF™ are trademarks or registered trademarks of Vitalograph Ltd.



Print Ref: 56027_1