



Matrix Park, 1 Western Avenue, Buckshaw Village, Chorley PR7 7NB
T : +44 (0) 1772 299 900 F : +44 (0) 1772 299 901

PRODUCT SPECIFICATION

Product Range		Conti® Patient Cleansing Wipes	
Product Name		Conti® Soft Large	
Product Code	Product Description	Size / Dimensions / Tolerances	Material / Composition / Active Ingredients
CSW110	Dry Wipe	320 x 280 mm (+/- 5%)	Material: 80% Polypropylene 20% Viscose Material: Weight: 30 GSM Wipe fold: C Fold Active ingredients: Not applicable
Intended Use		Patient dry wipe, suitable for everyday use.	
Classification		Unregulated	
Product/Device Single Use		Yes – individual wipes are single use.	
Components Latex Free		Yes	
pH Information		Not applicable	
Shelf Life		5-year shelf life	
Storage		Keep dry and protect from direct sunlight	
Sterilisation Method		Not applicable	
Special Instructions			
Do not flush, do not macerate.			
Pack Qty, Carton Qty, Case Qty, Pack/Pouch content Weight, Case Weight		Pallet Dimensions	
Pack: 100 (+/- 2%) Case: 32 Case Weight: 7 KG approx.		1200 X 1000 X 1380 mm	
Pallet Type	Pallet Configuration	Pallet Weight	
4-way	Cases per layer: 5 Layers per Pallet: 3 Quantity per pallet: 15	125 KG approx.	
Packaging Information and Size			
Primary packaging	Secondary Packaging	Tertiary Packaging	
PE Flow-wrap Film 40 micron thickness, Straight Perforation	Cardboard case (C5) Case dimensions: 626 x 390 x 418 mm	Not applicable	
Labelling Specification			
Primary (inner):			
Brand logo, product description, disposal icons, number of wipes, business logo, made in Britain logo, product code, business address and contact details, individual details box, artwork reference number, barcode, lot number.			
Secondary (middle):			
Business logo, business address and contact details, product code, case qty, lot number, keep away from sunlight symbol, keep dry symbol, barcode.			
Tertiary (outer):			
Not applicable			
Standards of Compliance			
Not applicable			

It is the responsibility of the user to verify that are using the current released version
HARD COPY UNCONTROLLED UNLESS STAMPED, SIGNED AND DATED BY AUTHORISED PERSONEL



Matrix Park, 1 Western Avenue, Buckshaw Village, Chorley PR7 7NB
T : +44 (0) 1772 299 900 F : +44 (0) 1772 299 901

Approval – Vernacare Limited	
Marketing	
Name: Elizabeth Jamieson	Signed: 
Position: Product Manager	Date: 16 Jan 2019
Quality	
Name:  R. Downham 16.1.19	Signed: 
Position: PRODUCT MARKETING DIR	Date: 16 JAN 19