

# Digital flat scale with Bluetooth interface and high capacity

**seca**  
Precision for health

## seca 813 bt

- + Solid steel frame, up to 200 kg capacity
- + Non-slip rubber base for safe standing during measuring
- + Large LCD digits for easy reading
- + Transmits weight via Bluetooth
- + Compatible with Google Fit and Apple Health for simple weight trend tracking



seca 813 bt

## Digital flat scale with Bluetooth interface and high capacity



### Pleasant and comfortable weighing

All platforms are not created equal. For this weighing platform seca has chosen a high quality, two-component raised profile material which is soft like rubber and pleasant to the touch. But don't be fooled – this material is also tough. It fulfills two requirements at the same time: comfortable weighing and an easy-to-clean surface that stands up to all commercially available cleaning and disinfecting products.

### A robust, sturdy scale – with a touch of elegance

seca has designed the seca 813 bt with a Bluetooth interface to transfer the Auto-HOLD value automatically to an external device. Thanks to its solid steel frame, the seca 813 bt can withstand high weight and features large LCD digits for easy reading. Furthermore, its wide platform with a low-profile facilitates trouble-free weighing.



### Paying attention to details

Weight and girth are no problem for the seca 813 bt. An issue with many scales is the display. Thanks to the 29 mm high LCD digits, it is extremely easy to read the measurements. Additionally, weight values can be shown live on the connected Bluetooth app if viewing the scale display is troublesome. You see, a scale from seca makes life more simple by paying attention to even the slightest details.

Capacity	200 kg
Graduation	100 g
Dimensions (W x H x D)	433 x 47 x 373 mm
Net weight	2.9 kg
Power supply	Batteries
Interfaces	Bluetooth
Functions	Tap-on automatic switch-on, Auto-HOLD, automatic switch-off, kg/lbs switch-over, Bluetooth connectivity
Scope of delivery	4 x AAA batteries

We reserve the right to make modifications.

seca united kingdom  
40 Barn Street  
Birmingham B5 5QB  
United Kingdom

T +44 121 643 9349  
F +44 121 633 3403  
E info.uk@seca.com

Precision for health  
seca.com

201181305021