



GIMA

MD SCOPE VIDEO OTOSCOPE - ELITE PACK

Code: 32160
Category: Video otoscopes/dermatoscopes
Unit of sale: 1 pc.
Minimum order: 1
Type: Medical device
Class: I



EAN13: 8023279321609

Description: MD Scope with:
- 1 camera probe 47 mm Ø 4.2 mm
- RCA video cable
- 24 disposable specula
- 2 AA alkaline battery
- carrying case 25.5 x 16 x 7 cm

Instant View

A clear and sharp image can be acquired and displayed momentarily on a builtin screen. Just press a single button and wait for 3 seconds.

Real Time Interaction

The medical image can be captured and shared simultaneously with patients. A true improved interaction between the health professionals and the patients.

Electronic medical record

A medical image is automatically generated to support most of the stand-alone health information system. Simply plug the cable into the video-out socket to retrieve the image.

Detachable camera modules are designed for different examinations as needed.

A video port provides simultaneous image to display on TV/monitor or to download to Mac/PC system.

- Image freeze function
- 2.4" colour LCD display
- Video out port
- Image resolution: 320 x 240 pixel
- Lightweight: 218 g
- Battery lifespan: 2 hours
- Size: 7.5 x 6.3 x 21.2 cm