



RELIABLE PROTECTION FOR YOUR AED



Dear customer,

This brochure shows our ARKY product collection. At ARKY, we believe in great quality, competitive prices and providing an excellent service. AED accessories are just as important as the AED device itself. For example the cabinets, wall mounts and backpacks.

By supplying you with these products, we ensure that AED's can be installed at any place and in any circumstance. In that way we increase the findability and mobility of AED's.

That is our philosophy and it drives us to improve our level of services day after day.

Since the beginning in 2008, we have proudly delivered our products to tens of thousands of people in over 30 countries worldwide. We are excited to see what the future holds, as we have only just got started!

Yours sincerely,

Rogier Balvers
Business Unit Manager at ARKY





ARKY POLYSTYRENE OUTDOOR CABINET



High quality

- Made of High Impact Polystyrene (HIPS)
- ▶ Protects the AED against all types of weather
- ▶ 100% dust- and waterproof
- Well sealed with brushed stainless steel hinges and closures
- Including mounting materials



Dimensions

- ▶ Inside dimensions: 28 x 32 x 18 cm
- Outside dimensions: 34 x 40 x 20 cm



Options

- Acoustic alarm
- Company branding
- Including breakaway seal
- Two heating elements (24V)
- Padlock



ARKY outdoor cabinet horizontal



ARKY outdoor cabinet vertical

Product no.	Description	Aluminum inlay	Alarm	Lockable (With breakaway seal/padlock)	Heating elements (24V)
60.212	ARKY outdoor cabinet without heating elements	•	•	•	
60.213	ARKY outdoor cabinet with heating elements	•	•	©	©



ARKY METAL INDOOR CABINET



High quality

- > Strong metal cabinet creates a seamless look
- ▶ Ideal for office buildings and public spaces
- Perfectly powder coated
- ▶ Logo engraved aluminum hand grip
- Including mounting materials



Dimensions

- ▶ 38 x 38 x 20 cm
- Fits defibrillators of all brands



Options

- Acoustic alarm
- Available in green and white
- ▶ ILCOR AED logos



ARKY white and green indoor cabinet



ARKY green indoor cabinet with ILCOR AED logo

Product no.	Description	Alarm	Fits defibrillators of all brands
60.111	ARKY white indoor cabinet		©
60.112	ARKY white indoor cabinet with alarm	•	©
60.121	ARKY green indoor cabinet		
60.122	ARKY green indoor cabinet with alarm	•	



ARKY PERSPEX WALL MOUNT



High quality

- Made of high quality Perspex
- Easy to mount
- Polymer length optimized for both transparency and sturdiness
- Mounting materials included



Dimensions

- ▶ 24 x 46 x 15 cm
- Specially designed for compact AED's



Options

Customized AED stickers



ARKY Perspex wall mount



ARKY Perspex wall mount - side

Product no.	Description	Including stickers	Including mounting materials
60.511	ARKY Perspex wall mount	•	•



ARKY AED BACKPACK



High quality

- Designed in collaboration with First Aid Services
- Weather and water resistant
- Robust character
- High visibility



Dimensions

- Compact model: 31 x 40 x 15 cm
- ▶ Large model: 38 x 56 x 19 cm



Options

- Both models are able to equip defibrillators of all brands
- Fast response time in case of an emergency
- Flexible storing solutions



ARKY AED backpack – backside



ARKY Compact and Large AED backpack

Product no.	Description	Weather and water resistant	Fits defibrillators of all brands
60.611	ARKY Compact AED backpack	•	•
60.612	ARKY Large AED backpack	•	©



ARKY SAFESET



Content

- Microfiber tissue
- Rescue blanket
- Pocket mask
- Scissors
- Razor
- Gloves




Purpose

▶ Essential accessories for CPR



Content of ARKY Safeset

Product no.	Description
70.112	Safeset



OUR NETWORK

Since the start in 2008 we delivered our products to over 30 countries worldwide. And we are constantly looking to further expand our business network. In case you are interested in our products, please do not hesitate to contact us. We will be pleased to inform you.

Contact

ARKY
Verlengde Bremenweg 12
9723 JV Groningen
The Netherlands

+31 (0)50 31 16 085

sales@arkycase.com

www.arkycase.com

