

STAINLESS STEEL GONIOMETERS



• **24900 FINGER GONIOMETER 6" - 15 cm**
Finger goniometer measures range of motion of finger joints. Linear scale in inches and centimeters on arms. 180° scales marked in 5° increments.



• **24901 180° DIGIT GONIOMETER 6" - 15 cm**
Allows taking readings using only one hand. Measures 110° flexion through 40° hyperextension in 5° increments.



• **24903 360° GONIOMETER 14" - 35 cm**
This goniometer has two opposing 180° scale and one 360° scale all marked in 1° increments. Sturdy adjustable knob, controlling tension of arms.




• **24902 180° ROBINSON GONIOMETER 7" - 17.5 cm**
Pocket size 180° Robinson goniometer for doctor and therapist. Linear scale marks 5° increments.




• **24904 180° EXTRA LONG ARMS GONIOMETER 14" - 35 cm**
Extra long arms, useful for measuring all joints. Arm has a linear scale marked in 1° increments.

GONIOMETERS


Used by orthopaedic surgeons to measure the range of movement of joints: elbows, knees and fingers.



• **27340 GONIOMETER - 205x45 mm**
With circular scale used to check the size of hip joints and the angle of the femoral neck on X-rays. Personalisation with your logo for minimum 2,500 pcs.



• **27339 GONIOMETER WITH PAIN SCALE RULER - 203x45x7 mm**



• **27338 GONIOMETER WITH ARMS (tolerance 1°)**
Measures analytical motility of various body segments in degrees (from 0 to 360 clockwise or anti-clockwise direction). Consists of two perspex arms (3 mm thick) that will remain open at any given position due to the friction at the joint between them. Two graded scales: one increasing scale and one decreasing scale.

BMI TAPE



• **27341 BMI - Tape Measure**
Pocket retractable tape measure 150 cm. BMI (Body Mass Index) calculator on front panel.

• **27342 BMI - Tape Measure**
Economical version of 27341.



27341

BODY TAPE




• **27343 BODY TAPE MEASURE**
Track your body measurements and progress. Measures any body part: arm, thigh, calf, chest, waist, hip... Measures in inches and centimeters; total measurement length: 60 inches / 152 cm. Sturdy vinyl tape.

Push-button retraction and locking feature ensures easy measurement and superior accuracy.



BMI CALCULATOR



• **27336 BMI CALCULATOR - box of 10**
Used in hospital and surgeries gives an indication of general health through the ratio of weight to height. Made of 2 rotating discs: 108 mm diameter (base) and 76 mm (top). Can be personalized with your logo for a minimum of 1,000 pieces.

PROFESSIONAL BODY - HEAD TAPE



- **27347 PROFESSIONAL BODY TAPE**
- **27350 PROFESSIONAL BODY TAPE with pediatric gauge box**

Patented, professional body tape measure approved by the Italian Ministry of Health for adult and pediatric use. Tape is coloured on both sides to help to identify metabolic risk in abdominal measure for male and female. Length of tape: 8-150 cm. Thanks to its rounded shape, it can be used also as a head measure tape. Packed in single box.

Body tape with pediatric gauge box (27350)
Specific tape for pediatric use. Supplied in a box with a ruler to compare the measured head circumference with the standard percentiles of children's growth from 0 to 36 months. Suitable for both sexes.

Level of metabolic risk for male (female on back)



27350

Multilanguage manual: GB, FR, IT, DE, ES, GR, PT, Arabic.