

## SECA MECHANICAL SCALES

## • 27241 SECA 750 SCALE

Mechanical flat scale with sturdy, powder coated steel casing. Pretty and stable with its white powder coated lacquer.

A hard wearing synthetic white leather covers weighing platform.

User manual in: GB, FR, IT, ES, DE, PT, PL.

**seca**

## • 27280 SECA 761 SCALE - medical use - Class III

## • 27285 SECA 762 SCALE

Traditional model with large easy-to-read dial.

Pointer professional scale in strong metal structure white painted. Consistent accuracy on any surface. Spacious and stable platform with no slip rubber mat. Tare function.

User manual in: GB, FR, IT, ES, DE, PT, PL.

27280

## TECHNICAL SPECIFICATIONS

Capacity: 150 kg - 330 lbs  
Sensitivity: 500 g  $\pm$  2 g (27285) - 1,000 g (27280)  
Size: 305x120x470 mm Weight: 3.5 kg

MEDICAL  
USE

27241

## TECHNICAL SPECIFICATIONS

Capacity: 150 kg - 330 lbs  
Size: 317x117x477 mm  
Sensitivity: 1 kg  
Weight: 3.5 kg

## SECA DIGITAL SCALES

## • 27265 SECA 813 DIGITAL FLAT SCALE

Electronic scale with high capacity up to 200 kg. Extra-wide platform and flat shape facilitate weighing of heavyweight patients and make stepping on very easy. Made of solid steel frame with tough rubber coating and non-slip rubber base, ensuring safety and stability. Kg/lbs/sts switch. Auto-hold, step-off and automatic switch ON/OFF functions. Multilingual manual: GB, FR, IT, ES, PT, DE, PL, NL, GR, DK, SE, FI, NO, TR, JP, RU.



27265

UP TO 200 KG

**seca**

## TECHNICAL SPECIFICATIONS

Capacity: 200 kg - 440 lbs  
Size: 433x373x47 mm  
Weight: 2,9 kg  
Sensitivity: 100 g - 0,2 lbs  
Platform size: 400x300x45 mm  
Batteries: 4xAAA 1,5V (not included)

## • 27257 SECA 807 DIGITAL SCALE

Digital scale with advanced technology for enhanced comfort. Extra-flat design to easily step onto the platform. Four specially constructed load cells made of stainless steel, guarantee high precision. Supplied with power sensors, extra-large LCD numbers and high-quality glass surface. Safety shock resistant glass construction. Step-off, auto-hold and automatic switch ON/OFF functions.



CE

27257

Multilingual manual: GB.

FR, IT, ES, PT, DE, PL, NL, GR, DK, SE, FI, NO, TR, JP, RU.

## TECHNICAL SPECIFICATIONS

Capacity: 150 kg  
Size: 310x310x37 mm  
Weight: 2 kg  
Sensitivity: 100 g  
Batteries: 4xAAA 1,5V (not included)

## iHEALTH FIT BODY ANALYSIS SCALE

See the list of compatible devices on [ihealthlabs.eu/support](http://ihealthlabs.eu/support)

FDA CE



iHealth  
MyVitals

23503

## • 23503 iHEALTH FIT HS2S WIRELESS BODY ANALYSIS SCALE

Weight and body composition monitoring on iPhone, iPad and Android smartphone.

- light and compact body analysis scale
- measures weight, BMI (Body Mass Index), lean, bone, muscle, water mass, body fat, daily caloric intake and visceral fat rating
- luminous LED display - 8 users
- track data, graphs and history on the iHealth MyVitals App
- store up to 500 measurements per user before synchronization
- share your medical data with your healthcare professional Bluetooth Sync, easy to install. Secure and free data storage.

## TECHNICAL SPECIFICATIONS

Connectivity: Bluetooth 4.0  
Dimensions - weight: 306x306x26.6 mm - 1,850 g  
Battery: 4 AAA batteries 1.5 V (included)  
Range - sensitivity: - weight 5-180 kg (11-400 lb) -  $\pm$ 1%  
- body fat 5%-65%  $\pm$ 1%

## GIMAFIT BLUETOOTH SCALE

## • 27086 GIMAFIT BODY SCALE with Bluetooth 5.0

Multifunction scale, for measuring body fat, water, muscle and bone rate and mass (BMI, BMR). Equipped with 4 precision G sensors. Platform material: glass + stainless steel. Domestic use.

**GIMA**

Bluetooth CE

## TECHNICAL SPECIFICATIONS

Capacity: 180 kg / 396 lb - minimum 0.15 kg  
Sensitivity: 200 g / 0.4 lb (>50 kg <100 kg)  
Size: 300x300x23 mm  
LED display: 80x40 mm white light  
Batteries 4xAAA (included)  
Switchable kg/lb

Multilingual manual in: GB, FR, IT, ES, PT, DE, PL, GR.



0.15-180 kg

27086