

SPECIFICATION

Cleartest Drug Test Hair Drug hair test



The CLEARTEST drug test hair is a chromatographic immunoassay for qualitative detection of 4 drugs (Cotinine, Methamphetamine, Morphine, Cocaine) and their metabolites in human hair.

The advantage of the drug hair test is the detection over a long period of time, which is not possible with any other sample material.

Evaluation after 5 minutes

Storage at 2 - 20 °C, do not freeze

Quick Info



HAIR



AFTER 5 MIN

SPECIFICATION

Cleartest Drug Test Hair

Type	PZN	REF	Unit
Drug hair test		C3 10400-1	1 test
Drug hair test		C3 10400-5	5 tests

Product	REF	Quantity	Unit price	Total

_____ plus tax

Personal details:

Customer number: _____

Company: _____

First name: Last name: _____

E-mail: _____

Street, street#: _____

Postal code, city: _____

I hereby place a binding order:

Company stamp / Address

Date / Sign

You are ordering under our terms and conditions, our valid terms and conditions can be found at <https://www.servoprax.de/en/terms-and-conditions/>