

## SPECIFICATION

### Allewyn Life Smith & Nephew



ALLEVYN Life dressing is a multi-layered design incorporating hydrocellular foam, hyper-absorber lock away core and masking layer which has been designed for people and their everyday life through.

Unique quadrilobe shape with a wide border designed to fit the contours of the human body more securely and allows patients to shower

The protective masking layer minimises the visual impact of strikethrough (in-vitro), providing the patient with confidence that their dressing may not attract the negative attention of others

Excellent fluid locking under pressure.

## SPECIFICATION

### Allevyn Life Smith & Nephew

Dimensions	OC	PZN	REF	Unit
10,3 x 10,3 cm	100 pcs.	09634019	L6 66801067	10 pcs.
12,9 x 12,9 cm	50 pcs.	09634031	L6 66801068	10 pcs.
15,4 x 15,4 cm	60 pcs.	09634048	L6 66801069	10 pcs.
21 x 21 cm	40 pcs.	09634060	L6 66801070	10 pcs.

Product	REF	Quantity	Unit price	Total

\_\_\_\_\_ plus tax

**Personal details:**

Customer number: \_\_\_\_\_

Company: \_\_\_\_\_

First name: Last name: \_\_\_\_\_

E-mail: \_\_\_\_\_

Street, street#: \_\_\_\_\_

Postal code, city: \_\_\_\_\_

**I hereby place a binding order:**

Company stamp / Address

\_\_\_\_\_  
Date / Sign

You are ordering under our terms and conditions, our valid terms and conditions can be found at <https://www.servoprax.de/en/terms-and-conditions/>