

ECOLOGICAL THERMOMETERS - MERCURY FREE



• 25551 ECOLOGICAL OVAL THERMOMETER - box of 12
 • 25587 SHAKER for ecological thermometers - box of 24
Mercury free clinical thermometer with gallium alloy instead of mercury. Safe, precise and functional. Glass structure. Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, RO, SK, GR, Arabic.

• 25586 ECOLOGICAL OVAL THERMOMETER with shaker - box of 1
 • 25589 ECOLOGICAL OVAL THERMOMETER with shaker - display box of 24 pcs.
 Same as 25551 but supplied with a shaker to help shake down gallium after use.

WIRELESS THERMOMETER BLUETOOTH 4.0

• 25548 TEMP SITTER - baby wireless thermometer
 Accurate and innovative thin thermometer. Transmission to Temp Sitter App (iOS, Android). To be placed under baby armpit. Update frequency every 3 seconds. High temp alarm can be set in the App. Waterproof and material of safe quality. Supplied with 10 adhesives to fix thermometer to baby armpit, manual (GB, FR, IT, ES).



TECHNICAL SPECIFICATIONS
 Battery life: 3,000 hours, not replaceable
 Measurement range: 25-45 °C (±0.05 between 35-38.5 °C)
 Size: 51x31.8x8.5 mm - Weight: 9.6 g - Switchable °C/°F

OMRON DIGITAL THERMOMETERS

- 25603 ECO TEMP BASIC THERMOMETER - MC 246-E
- 25604 FLEX TEMP SMART THERMOMETER - MC 343F-E



	TECHNICAL SPECIFICATIONS	
Model	25603 (Eco Basic)	25604 (Flex Smart)
Measurement Type	Oral, Armpit, Rectal	Oral, Armpit, Rectal
Flexible tip		•
Beeper (signals when a measurement is complete)	•	•
Memory	last reading	last reading
°C/°F selection	•	•
Large, 3 digit display	•	•
Water resistant	•	•
Automatic shut off	•	•

LIQUID CRYSTAL FOREHEAD THERMOMETERS

- 25599 FEVERTEST - bag of 10 blister
- 25596 FEVERTEST BULK - box of 100 - hospital use
- 25600 FEVERCONTROL - bag of 10 blister



25599

25600

48h FEVERMONITOR

- 25601 48h FEVER MONITOR - bag of 10 blister of 3 pcs.



25601

Innovative "continual" fever indicator for 48h monitoring, for babies. Disposable, non toxic, mercury free, safe, unbreakable, reliable, easy to read, liquid crystal thermometers. It indicates temperature in 30 seconds. Medical Device.