

FIRST AID PLASTIC CASES AND METAL CABINETS



34196



34197

- 34196 MEDIUM KIT - plastic case Size: 275x240x123 mm - 2 kg
- 34197 MEDIUM KIT - metal cabinet Size: 370x300x140 mm - 3.8 kg
- 34187 REFILL PACK for code 34196-97
- 34141 EMPTY METAL CABINET 370x300x140 mm



34198



34199



34195

- 34195 LARGE KIT - plastic case 1 Size: 460x345x135 mm - 6.2 kg
- 34198 LARGE KIT - plastic case 2 Size: 395x270x135 mm - 6.2 kg
- 34199 LARGE KIT - metal cabinet Size: 460x300x140 mm - 7.5 kg
- 34189 REFILL PACK for code 34195-34198-34199 without sphygmomanometer
- 34191 EMPTY PLASTIC CASE 1 460x345x135 mm
- 34138 EMPTY PLASTIC CASE 2 395x270x135 mm



- 34143 METAL CABINET
Supplied with 3 internal metal shelves, removable bottle floor support, 2 small shelves on the door, lock and key. Size 402x582x202 mm - 12.9 kg
- 34145 EMPTY METAL CABINET 402x582x202 mm

STANDARD CASES / CABINETS CONTENT	34194 34184-5	34196 34197	34195 34198-9	34143	34151	SPARES see page
Bags liquid soap		3*		3		418
Bags antiseptic wipes		3*		3		418-419
Bottle antiseptic solution 125 ml	1	1				-
Bottle antiseptic solution 250 ml				1		419
Bottle antiseptic solution 500 ml			2	2		-
Oxygenated water bottle				1	2	419
Ammonia bag		2*		2		-
Pair protective gloves, sterile	2	2	5	5		28-29
Box of 10 assorted plasters	1	1	2		5	410
Box of 100 assorted plasters Plastosan				1	2	410
Box of 100 7x2 cm Plastosan				1		410
Adhesive plaster roll 2,5 cm x 5 m	1	1	2	2		411
Elastic bandage with clip 4,5 cm x 7 m				1	2	407
Pack of 8 assorted gauze bandages				1		407
Gauze bandage 10 cm x 3,5 m	1	1				407
Tubular net gauze			1	1		403
Pack of sterile compress 10 x 10 cm	3	3	10	10	1	401
Pack of sterile compress 18 x 40 cm	1	3	2	6	10	401
Pack of cotton wool	1	1	1	3	2	403
Mask with face shield			1	1	1	17
Pair of scissors cm 14,5	1	1	1	1	1	374
Sterile forceps	1	1	2	2	1	50
Haemostatic lacet	1	1	3	3	2	76
Sphygmomanometer + Stethoscope			1	1	3	275-310
Sterile blanket 40x60 cm		1	2	2	1	-
Pack of instant ice, single use	1	1	2	2	2	380
Digital thermometer			1	1	2	102
Trash bag	1	1	2	5	1	-
Irrigating solution 250 ml	1	1			1	414
Irrigating solution 500 ml			3	3		-
Triangular blanket 96x96x136 cm		1		4	3	-
Isothermic blanket gold/silver 160x210 cm				1		380
Sterile burn gel sachet 3.5 g		1*		3		380
Compressive Emocontrol bandage				1		407
Pinocchio helicoidal open mouth Vento				1		370
Autodavable resuscitator with facemask N° 4						371
Oxygen bottle 0.5 l empty - UNI standard					1	369
Guedel airways					1	372
Atoxic tube for oxygen					3	374
Mouth opener					1	370

*Only 34197

• 34193 MULTIRED LARGE KIT

- plastic case 3
size: 443x338x147 mm
Made of strong triopolymer plastic which is shiny smooth and easy to clean. Neopren seal fitted to make it dustproof. A wall bracket secures it on any house wall or booth. Internal separators and safe closure with 2 swivel catches. Weight: 3.3 (empty), 8.4 kg (filled).



34193

STANDARD ACCESSORIES for 34193

5 pairs of sterile gloves	3 tourniquets
1 safety mask with eye protection visor	2 disposable instant ICE packs
3 sterile saline solutions - 500 ml	2 sanitary bags
2 Iodopovidone disinfectant solutions 10% - 500 ml	1 clinical digital thermometer
10 sterile gauzes in sachets 10x10 cm	1 sphygmomanometer with stethoscope
2 sterile gauzes in sachets 18x40 cm	2 gauze bandages 10 cm x 3.5 m
2 sterile bandages 40x60 cm	1 isothermic blanket gold/silver 160x210 cm
1 pack of cotton	1 triangular TNT bandage 136x96x96 cm
2 sterile tweezers	1 pack of 8 sachet soap/disinfectant
1 pack of net elastic bandage	2 elastic bandages 4 m x 6 cm
2 adhesive plasters 5 m x 2.5 cm	1 sterile compressive bandage 80x100
2 Plastosan 20 assorted dressings	1 Pinocchio mouth opener with airways
1 pair of lister scissors 14.5 cm	2 burn gel sachets 3.5 g