


## SURGICAL SKIN MARKER - STERILE


• **33181 P3 SKIN MARKER® - fine tip - sterile - box of 10**

The P3® Surgical Skin Marker features a versatile tip for precise fine line marking preferred in plastic surgery or when held at an angle, broad line marking for applications in general, neuro and vascular surgery, as well as radiology and radiation therapy. The marker's high visibility Gentian violet ink is non-toxic, non tattooing and surgical prep-resistant. Markings can be removed from the skin with isopropyl (rubbing) alcohol. The package is easy to open, HIBC bar coded and includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Sterile EO. Individually packed.



33181

**TWO TIPS IN ONE:**




For heavier lines, hold the pen at an angle and use the shoulder of the tips (2 mm)



Hold the pen vertically while applying a fine line (0.2 mm)



## SURGICAL SKIN MARKERS - STERILE




33182

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, DK, PL, CZ, FI, GR.

• **33182 SURGICAL SKIN MARKER - single tip - sterile - box of 100**  
Single tip marker ideal for plastic surgery. Hypoallergenic and violet non-toxic ink. Marking can be removed by alcohol and liquid soap. Intended for single use. Includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Individually packed.

TECHNICAL SPECIFICATIONS	
Tip size:	1.0 mm
Diameter:	9.5 mm
Length:	167 mm
Ink colour:	violet
Washable:	no



33183

• **33183 SURGICAL SKIN MARKER - double tip - sterile - box of 100**  
Double tip marker ideal for plastic surgery. Hypoallergenic and violet non-toxic ink. Marking can be removed by alcohol and liquid soap. Intended for single use. Includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Individually packed.

TECHNICAL SPECIFICATIONS	
Tip size:	fine tip: 0.5mm
Diameter:	11 mm
Length:	140 mm
Ink colour:	violet
Washable:	no

## SKIN MARKER

• **33179 SKIN MARKER - dual tips gentian violet - box of 10**  
Double tip skin marker with polypropylene cap and body. Polyester fiber tips.  
**Ultra fine tip utility marker** to write on adhesive labels and infusion bags.

Multilingual box and manual: GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arabic.



33179

**Regular tip surgical skin marker** for use on patient skin. Washable with water and soap. Size 138xØ 12.5 mm.



With ruler in mm

33179

Ultra fine tip

Regular tip

## DERMATOGRAPH MARKERS

- **33175 DERMATOGRAPH MARKER - blue - box of 6**
- **33176 DERMATOGRAPH MARKER - white - box of 6**
- **33177 DERMATOGRAPH MARKER - red - box of 6**
- **33178 DERMATOGRAPH MARKER - mixed colours - box of 6\***

\*2 pencils each colour

Special dermatograph pencil to mark the skin for identification of areas to be treated.

Made of cedar wood with wax lead. Size 175xØ 7 mm. Suitable for general surgery, dermatology, orthopaedy, radiology, ...

Pencils are sold non sharpened in box of 6 with 1 pencil sharpener included.

Multilingual box and manual:

GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arabic.



33175

33176



33177

33178

