

## SPECIFICATION

### Omron HBP-1120-E



Professional upper arm blood pressure monitor for auscultation mode and oscillometric measurement with zero indicator.

- Compact design: fits on any desk
- Simple user interface: 3 buttons and large display
- Display: 7-segment LCD
- Compatible with 5 cuff sizes (from 12 to 50 cm), cuffs can be disinfected
- Detects irregular heartbeats and body movements
- Includes M cuff (22 - 32 cm), AC adapter (power supply unit) and instruction manual

Device dimension without cuff (WxDxH): 130 x 175 x 120 mm

## SPECIFICATION

### Omron HBP-1120-E

Type	PZN	REF	Unit
Omron HBP-1120-E	14057535	E3 HBP-1120E	1 piece

Product	REF	Quantity	Unit price	Total

\_\_\_\_\_ plus tax

**Personal details:**

Customer number: \_\_\_\_\_

Company: \_\_\_\_\_

First name: Last name: \_\_\_\_\_

E-mail: \_\_\_\_\_

Street, street#: \_\_\_\_\_

Postal code, city: \_\_\_\_\_

**I hereby place a binding order:**

Company stamp / Address

\_\_\_\_\_  
Date / Sign

You are ordering under our terms and conditions, our valid terms and conditions can be found at <https://www.servoprax.de/en/terms-and-conditions/>