

PREGNANCY TESTS: MIDSTREAM, STRIP AND CASSETTE

SELF TESTS CE 0123

MIDSTREAM

29098

CE 0197

Test di gravidanza
 • Test de grossesse
 • Test de embarazo

Extra early detection 3 days before you miss your period

SHELF LIFE 36 MONTHS

CE 0123

29097

29103

Large wipe allows a better absorption of the urine sample reducing the risk of false results

CUSTOMIZED PREGNANCY SELF TESTS
 minimum 5,000 pcs.

GIMA offers a complete line of pregnancy tests, both for OTC market (self test) and for professional use. All products are packed in cardboard boxes with multilingual user manual (GB, IT, ES, FR, GR, PT, DE*, PL*). *Manual only

Main features:

- each test is individually packed in sealed aluminium pouch
- very quick, in 3 minutes
- early detection: it can be used since the first day you miss your period (three days before for 29103)
- accuracy: 99.9% (29090-2), 99.0% (29097-8, 29103)
- sensitivity: 20 mIU/ml (29090, 29091, 29092, 29097, 29098) 10 mIU/ml (29103)
- clinically tested

PROFESSIONAL TESTS CE

MIDSTREAM

29090

STRIP

29091

CASSETTE

29092



Read results in 3 minutes (29090, 29097, 29098, 29103)

Step 1: OR

Step 2: Until complete colour change of the indicator

Step 3: 3 minutes

Positive: T C T C

Negative: T C

Invalid: T C T C

PREGNANCY TESTS			
GIMA code	SELF TESTS	Type	Package
29097	SELF TEST 20 mIU/ml	midstream	box of 1 pc
29098	SELF TEST 20 mIU/ml	midstream	box of 2 pcs
29103	ULTRA SELF TEST 10 mIU/ml	midstream (large wipe)	box of 1 pc
PROFESSIONAL TESTS			
29090	PROFESSIONAL TEST	midstream	box of 25 pcs
29091	PROFESSIONAL TEST	strip 4 mm	box of 50 pcs
29092	PROFESSIONAL TEST	cassette	box of 40 pcs

ULTRA SENSITIVE OVULATION TESTS

OVULATION TESTS:
 Easy to use self test for determination of the highest peak of fertility in woman. Early detection test, 99.9% accurate and clinically tested. Sensitivity: 35 mIU/ml (29124), 30 mIU/ml (29125). All products are packed in single aluminium pouches in a cardboard boxes with multilingual user manual (GB, IT, ES, FR, GR, PT, DE, PL). Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months. Results in 3 minutes (29124) or 5 minutes (29125).

GIMA code	OVULATION TESTS	Type	Package
29124	SELF TEST	midstream	box of 5 pcs
29125	SELF TEST	strip 4 mm	box of 25 pcs

SELF TEST CE 0123

29124

29125

29124

29125

BASAL THERMOMETER


• 25608 DIGITAL BASAL THERMOMETER °C - hang box - water-proof
 With 4 digits, 2 decimals, for fertility control.

4 DIGITS, 2 DECIMALS CE 0476

25608

PREGNANCY DISC

• 29099 PREGNANCY DISC
 Forecasts the forthcoming date of birth. Used by obstetrician and practitioners. Useful for assessment of drug administration. 2 rotating discs riveted in the center: 108 mm diameter (base), 76 mm diameter (top). Can be personalized with your logo for minimum 1,000 pcs.



YOUR LOGO HERE
 minimum 1,000 pcs.