

BURNFREE® - BURN CARE PRODUCTS

A complete range of burn care products with superior benefits:

- dissipates heat on contact
- cools the burn site, relieves pain, calms the burn victim, helps stop the burn process
- covers the burn, helps to prevent contamination
- accommodates swelling or the need to flex joints thanks to elastic texture

Sterile burn and wounds dressings

Made of a soft open-cell, non-adherent, inert foam, saturated with Burnfree® gel, does not leave residues imbedded in the wound. Without lidocaine, which makes it suitable for children.

GIMA code	BURNFREE® PRODUCTS
34869	Burn Dressing - 5x5 cm - box of 120 pcs
34870	Burn Dressing - 5x15 cm - box of 10 pcs
34871	Burn Dressing - 10x10 cm - box of 10 pcs
34872	Burn Dressing - 20x20 cm - box of 5 pcs
34873	Burn Dressing - 30x30 cm - box of 20 pcs
34875	Pain Relieving Gel - bottle 120 ml
34883	Pain Relieving Gel - 3.5 g sachets - box of 20 pcs
34882	Pain Relieving Gel - 3.5 g sachets - box of 1,000 pcs
34878	Burn Kit (2 sterile dressing 5x15 cm and 2 10x10 cm, 1 gel bottle 120 ml, 20 gel sachets 3.5g and 2 rollers gauzes)

Pain relieving gel


Ideal for all minor burns (including sunburns), cuts and abrasions, available in bottle or single use dose 3.5 g.



EMERGENCY BLANKET

• 34086 SICURTEK EMERGENCY BLANKET

All-purpose emergency blanket. Waterproof, anti-putridic, protects from cold and/or heat.



Made of polypropylene Fralom fibre coupled to metallized polyethylene fibre.

Blanket colour availability may change in the time

Size: 200x150 cm

BURNFREE® - FIRE BLANKET

BURNFREE FIRE/TRAUMA BLANKETS

Fire/trauma blanket for first aid treatment of burn injuries. It provides physical protection from heat, fire and smoke. It draws the heat out to prevent further tissue and nerve damage.

- it stops the burning process
 - it removes heat from burn site
 - it extinguishes small fires
 - medical grade ingredients.
- Blanket rapidly cools and soothes the burn. Supplied with bag.



Bags with strap for easy transport (34880-1 only)

GIMA code	Blanket Size	Weight
34879	91x76 cm	2.3 kg
34880	183x152 cm	6.9 kg
34881	244x152 cm	10.9 kg

ISOTHERMIC COVER

• 34085 EMERGENCY COVER

Suitable for protection against cold and heat. Silver side in keeps the warmth of the body, gold side out reflects sunbeams.



Size: 160x210 cm
Weight: 63 g

Radio transparent

Instructions: GB, FR, IT, ES, PT, DE, HU, RO, SE, FI, GR, HR, Arabic

FIRE RETARDANT BLANKET

• 34087 FIRE RETARDANT BLANKET

Made of glass fibre. Can be used to extinguish small fires. Size: 150x180 cm Standard EN1869 1997.



INSTANT ICE: SPRAY OR BAG IN TNT OR PE



Instant ice with multilingual bag: GB, FR, IT, ES.



34110 TNT



34111 PE



34112

- 34110 INSTANT ICE - TNT 14x18 cm bag (non woven fabric)
 - 34111 INSTANT ICE - PE 14x18 cm bag (polyethylene)
- minimum order: pack of 25 bags. (38x25xh 17 cm)

Icy-pack "instant ice" is recommended for cryotherapy. Excellent in case of traumas to soothe the pain and to prevent haematomas, swellings, etc.

Easy and quick to use anywhere and anytime.

Press the bag to break the inner bubble and shake to mix contents. The blend you obtain drops temperature of the pack down. It reaches -4°C up on activation, it reaches 0°C in 20 minutes (at room temperature of 20 °C). Weight 250 g.

• 34112 SPRAY INSTANT ICE - bottle 400 ml

Useful for injuries, scrapes, falls or sprains, in first aid or sport applications. Spray on the injured area at 20 cm distance. Content: gas (isobutane, butane, propane), liquid (methylal, isopropyl, myristate, menthol).