

## LAERDAL LITTLE FAMILY PACK

### • 56512 LAERDAL LITTLE FAMILY PACK QCPR - 136-01050

Little Family Pack is a convenient package of age-specific CPR training manikins (Little Anne QCPR, Little Junior QCPR and Little Baby QCPR) with a wheeled carrying case for ease of transport and storage. All are provided with QCPR feedback technology to help instructors to improve CPR training quality, efficiency and learner engagement. In compliance with current guidelines, they provide data about the rate and depth of chest compression, complete pressure relaxation, the time of CPR pause and ventilation, through the new QCPR Training Apps, skill guide available via RC cable.

It comes complete with: **Little Anne QCPR Manikin** (manikin, 2 faces, 2 airways), **Little Junior QCPR Manikin** (manikin, 2 faces, 2 airways), **Little Baby QCPR Manikin** (manikin, 6 airways\*), 18 manikin wipes, 3 training mats\*\* and a large bag on wheels.

\*Included with 56512, not included with 56509, 56510 \*\*Only 1 with 56509, 56510

- 56509 LAERDAL LITTLE BABY QCPR MANIKIN - 133-01050
- 56510 LAERDAL LITTLE ANNE QCPR MANIKIN - 123-01050



56512  
3 manikins  
with carrying  
case



Little Family QCPR  
LAERDAL

Little Baby QCPR

Little Junior  
QCPR

Little Anne  
QCPR

## LAERDAL RESUSCI ANNE

### • 56523 LAERDAL RESUSCI ANNE QCPR FULL BODY + TROLLEY - 171-01260

Resusci Anne QCPR is built with realistic anatomy and is fully guideline compliant. It allows multiple feedback options providing opportunity to focus on learner competency: the learner will experience appropriate chest resistance, see chest rise during ventilation and need to perform head tilt, chin lift to ensure an open airway. The QCPR Training and skill report App Android/iPhone technology enables healthcare providers to obtain competence to deliver quality CPR in accordance with the latest guidelines using Bluetooth devices. Supplied with full body with installed 45kg compression spring, two extra compression springs (30 kg and 60 kg), manikin wipes, 3 face skins in paper bag, 2 airway assembly (standard lung), trolley case, USB 2A AC power adapter, USB cable A-C for the connection with the PC instead of Bluetooth and for battery charging (Li-ion, 2 cells). The connection of the skill guide takes place only via an RC cable connected to the shoulder of the manikin.

User guide (GB, FR, DE, ES, BR, NL, NO, SE, DK, FI, JP, KO, RU, PL).



56523  
Full body  
manikin  
with trolley  
case

Skill guide connected to manikin shoulder

## PRACTI-MAN CPR MANIKIN



Multilingual manual:  
GB, FR\*, IT, ES\*, PT\*, DE\*  
for 56500 and 56501  
\*Only on www.gimaitaly.com

56500

Includes: 1 Practi-man,  
5 lungs, 2 valves,  
1 transport bag with mat.

- 56500 PRACTI-MAN ADVANCE CPR MANIKIN
  - 56501 4 PRACTI-MAN ADVANCE CPR MANIKINS
- CPR manikins for massive training events:
- adult and child CPR training with a single manikin
  - 3-position selector dial: adult, child and test
  - real-like system for breathing function
  - jaw thrust maneuver
  - AED training
  - 56504 PRACTI-MAN LUNGS - box of 5 - spare
  - 56505 PRACTI-MAN VALVES - box of 2 - spare

Jaw thrust  
maneuver



Includes:  
4 Practi-man,  
20 lungs,  
8 valves,  
1 transport bag  
with mat



56501

## CHEST-eR REANIMATION DEVICE

### • 56530 CHEST-eR REANIMATION DEVICE

Chest-eR is a new generation device which assists and improves cardiac massage practice. It can be used by rescuers, hospital staff, common people and educators.

It has been designed to reduce the incidence of internal injuries, providing protection during chest compressions.

A special 3 layer internal structure helps reducing dangerous stresses produced by incorrect application of force, redistributing excessive forces and guaranteeing hygiene.

The device includes luminous indicators that inform and show the rescuer on correct pressure and compression rate.

The device perfectly combines with manikins updated according to International Guidelines (i.e. Laerdal Medical, Brayden, .....)

and can be used with defibrillators. Features:


- intuitive and easy-to-use
- compact and lightweight, used with 1.5V battery
- proper resuscitation in complete safety
- ensures maximum effectiveness of cardiac massage
- non-Newtonian materials dissipates impact of too violent compressions absorbing up to 90% of energy

Medical Device Class I. Supplied with 2 protection sheaths and multilingual manual (GB, FR, IT, ES, DE)

### • 56531 PROTECTION SHEATH - box of 30

It guarantees absolute hygiene and prevents risk of infections.

Rigid middle layer allows  
forces distribution in the  
lower sternum area



56530

Luminous  
electronic  
feedback