

SPECIFICATION

mikrozid AF wipes Jumbo



Practical disinfectant wipes for rapid disinfection of medical equipment.

- Microbiological effectiveness: bactericidal, tuberculocidal, levurocidal, fungicidal, limited virucidal, norovirus, polyoma SV40, rotavirus, adenovirus (type 5).
- Virucidal in accordance with RKI-recommendation 01/2004
- IHO-listed, VAH certificate

Notes

Danger warnings

Signal word: Warning

Danger warnings

H226 - Flammable liquid and vapour.

H319 - Causes serious eye irritation.

H336 - May cause drowsiness or dizziness.

Safety Instructions

P102 - Keep out of reach of children.

P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280 - Wear protective gloves/eye protection.

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 - If eye irritation persists: Get medical advice/attention.

Additional warnings

This substance/mixture does not contain any components at concentrations of 0.1% or higher that are classified as either persistent, bioaccumulative and toxic (PBT) or very persistent and very bioaccumulative (vPvB). Environmental information: The substance/mixture does not contain any components that exhibit endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated Regulation (EU) 2017/2100 or Commission Delegated Regulation (EU) 2018/605 in quantities of 0.1% or more. Toxicological information: The substance/mixture does not contain any components that exhibit endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated Regulation (EU) 2017/2100 or Commission Delegated Regulation (EU) 2018/605 in quantities of 0.1% or more. Has a degreasing effect on the skin.

Biocide Registrationsnummer:N-40580

Use biocidal products with caution. Always read the label and product information before use.

?

SPECIFICATION

mikrozid AF wipes Jumbo

Type	OC	PZN	REF	Unit
Jumbo Dispenser with 220 towels 20 x 20 cm	10 pcs.	17386989	E2 70003093	1 can
Jumbo Refill with 220 towels 20 x 20 cm		17386995	E2 70003094	12 x 220 towels

Product	REF	Quantity	Unit price	Total

_____ plus tax

Personal details:

Customer number: _____

Company: _____

First name: Last name: _____

E-mail: _____

Street, street#: _____

Postal code, city: _____

I hereby place a binding order:

Company stamp / Address

Date / Sign

You are ordering under our terms and conditions, our valid terms and conditions can be found at <https://www.servoprax.de/en/terms-and-conditions/>