


# Visitec® Eye Shield

Single-Use Sterile Material



## DATA SHEET




<b>Product Class:</b> FDA Class I MDD Class I, sterile, single-use According to European Medical Device Directive 93/42/EEC	<b>Shelf Life:</b> <ul style="list-style-type: none"><li>• Five (5) years from manufacturing date.</li></ul>
<b>Manufacturer:</b> Beaver-Visitec International, Inc. Waltham, MA USA	<b>Intended use:</b> <ul style="list-style-type: none"><li>• Intended to be used to protect the eye pre or post surgery or with ocular trauma.</li></ul>
<b>Certification Information:</b> EC Certificate ISO 13485	<b>Product Composition:</b> <ul style="list-style-type: none"><li>• K-Resin, Lightweight Copolymer Plastic</li><li>• Absence of: Latex / PVC / Phthalates / Animal origin components</li></ul>
<b>For Sterile Product:</b> <ul style="list-style-type: none"><li>• Gamma Irradiation (cobalt 60)</li></ul>	<b>Packaging Material:</b> <b>Pouch:</b> <ul style="list-style-type: none"><li>• Polyethylene Film, Coated Tyvek Lid</li></ul> <b>Shipper:</b> <ul style="list-style-type: none"><li>• Paper Carton</li></ul>

**General Features:**

- Clear lightweight plastic. Features smooth edges and ventilation holes for patient comfort.

**Warnings and Precautions:** See box label and insert

**Product Information:**

<b>Visitec® Eye-Shield</b>		
<b>Item #</b>	<b>Description</b>	<b>Quantity</b>
<b>581019</b>	<p><b>Visitec® Left/Right Eye Shield 8.5 x 6.5cm (3 1/2 x 2 1/2 in)</b></p> 	25 Pairs
<b>581040</b>	<p><b>Visitec® Universal Eye Shield, 7.5 x 6.5cm (3 x 2 1/2 in)</b></p> 	50/box
<b>581708</b>	<p><b>Visitec® Universal Pediatric Eye Shield 5.1 x 4.1cm (2 x 1 5/8in)</b></p> 	25/box



**For US Customers Only:**

10 CityPoint  
 500 Totten Pond Rd  
 Waltham, MA 02451 USA  
 Office: 866.906.8080  
 Fax: 866.906.4304  
[www.beaver-visitec.com](http://www.beaver-visitec.com)

**All Other Customers:**

85c Park Drive, Milton Park  
 Abingdon, Oxfordshire  
 OX14 4RY UK  
 Office: +44 1865 601 256  
 Fax: +44 1865 595 761  
[www.beaver-visitec.com](http://www.beaver-visitec.com)