

CUBO EMERGENCY BAG



27163

• **27163 CUBO EMERGENCY BAG - PVC tarpaulin - red**
 Portable emergency bag with internal and external zip pocket and 2 removable pouches with transparent window. Interior made of strong, water resistant 600D polyester with outer PVC. Provided with double reflective yellow strip. Size: 28x34x13 cm.

LIFE-2 EMERGENCY BAG



27161

Shoulder or hand portable.

• **27161 LIFE-2 EMERGENCY BAG with 5 coloured pouches**
 2 section first-aid bag with double opening zip: bottom side with shockproof inner basin, upper side with 5 removable coloured pouches (model E1) with transparent window. Made of strong, water resistant 600D polyester PVC coated material with reflective yellow strips. Size: 49x28xh 33 cm.

EMERGENCY TROLLEY BAG

Trolley bag and rucksack made of strong, water resistant 600D polyester (27154) and PVC (27173) with durable cold resistant wheels.

• **27154 EMERGENCY SMART TROLLEY BAG - 56x32.5x49 cm**

Extremely capacious with internal adjustable compartments which permit personalized organization of the space. Provided with external pockets, internal divisions and one removable pouch with transparent window. Double yellow reflective strips. Waterproof rubber bottom to prevent water infiltration.



27154

• **27173 EMERGENCY LOGIC 2 PVC TROLLEY RUCKSACK - 49x31x57 cm**

Supplied with 6 coloured pouches with transparent window and velcro strips, 2 frontal pockets and 2 lateral pockets. Internal velcro strips and elastics, to retain instruments. Central division with zip and elastics. Provided with yellow reflective strips and shoulder straps.



27173



27154

27173

POUCHES WITH VELCRO

Padded back and front **E5 27197**



E3 27191



E2 27184



E4 27194



POUCHES WITH VELCRO MADE OF 600D POLYESTER

All pouches are provided with a large strip of female velcro on back or side to be fixed to Gima emergency bags.



E1 27180



E6 27200

| Model | Size - cm | RED | BLUE | YELLOW | GREEN | Standard with |
|-------|-----------|-------|-------------|--------|-------|----------------------|
| E1 | 20x11x8 | 27180 | 27181 | 27182 | 27183 | 27161 |
| E2 | 28x12x10 | 27184 | 27185 | 27186 | 27187 | 27170/1/4/5, 27150-8 |
| E3 | 28x15x11 | 27188 | 27189 | 27190 | 27191 | 27171 |
| E4 | 33x15x11 | 27192 | - | 27194 | - | 27151/2/7/8, 27175 |
| E5 | 29x19x11 | 27196 | 27197 | 27198 | 27199 | 27174 |
| E6 | 38x10x5 | 27200 | grey colour | | | |

Padded back

NYLON BAGS IN 5 COLOURS

| NYLON BAGS - 22x15.5 cm | | | | |
|-------------------------|-------|-------|-------|-------|
| Blue | Green | Red | Black | Grey |
| 32791 | 32793 | 32794 | 32795 | 32796 |



Provided with zip closure and external pocket with velcro closure