


# GIMA ECG THERMAL PAPERS FOR DIFFERENT ECG BRANDS



MADE IN ITALY

GIMA code	ECG MANUFACTURER AND MODEL	Type	Width mm	Length m	Sheet	Marker	Grid*	Min. order	GIMA code	ECG MANUFACTURER AND MODEL	Type	Width mm	Length m	Sheet	Marker	Grid*	Min. order
32951	ASPEL Ascard A4, B56, Mr. Blue, Mr. Silver 3 Topaz, grey, v.07.204, orange, red, v.07.303	roll	112	25	-	no	or/red	10	32966	LTD450, LTD455/S, LTD455/W	pack	210	295	150	yes	or/red	1
32996	Ascard B5 eco, Mr. Green, mint	roll	58	25	-	no	or/red	20	32963	SE-1	roll	50	25	-	no	or/red	20
32952	Ascard, Mr. Gold, Gold/3 BIO CARE	roll	210	25	-	no	or/red	5	32964	SE-3	roll	80	25	-	no	or/red	10
32968	101G	roll	50	20	-	no	or/red	20	32965	SE-601	pack	110	140	143	yes	or/red	20
32989	300G, 3010	roll	63	30	-	no	or/red	20	32966	SE-6, 12, 1200 EXP	pack	210	295	150	yes	or/red	1
32969	IE3, IE300	roll	80	20	-	no	or/red	10	32960	SE-6, 12, 1200 EXP	roll	210	30	-	no	or/red	5
33002	3030, 6010	roll	112	27	-	no	or/red	10	32962	SE/1201	pack	210	140	215	yes	or/red	10
32992	1200	roll	210	30	-	no	or/red	5	<b>ELETTRONICA TRENITINA/CARDIETTE/CARDIOLINE</b>								
32962	1210, 1215	pack	210	140	215	yes	or/red	10	32971	Autoruler 12/1, comp. 1x12	roll	50	30	-	yes	blue	20
33004	IE12A, IE-15, IE-20	pack	210	150	200	yes	or/red	10	32972	Start 100	roll	130	27	-	yes	blue	10
33005	ECG-1260	pack	210	140	140	orange	orange	10	32973	AR 600, AR 600 Adv	roll	60	15	-	yes	blue	25
33007	IE6, ECG 3030, ECG 6010	pack	112	140	142	orange	orange	20	32977	AR 600 Adv, AR 600	pack	60	75	250	yes	blue	25
32969	BLT, 1203A, E 30	roll	80	20	-	no	or/red	10	32974	AR 1200	roll	120	18	-	yes	blue	10
32992	BLT-1212	roll	210	30	-	orange	orange	5	32978	AR 1200 (View/Adv)	pack	120	100	300	yes	blue	10
32966	E70, E80	pack	210	295	150	yes	or/red	1	32975	Daedalus, AR2100	roll	210	20	-	yes	blue	5
32950	Cardiocare/Cardiotouch	roll	215	25	-	no	or/red	5	32976	Daedalus, Start 200, AR2100	pack	210	150	200	yes	blue	8
33046	BIOSET-HORMAN Bioset 3500	roll	110	30	-	orange	orange	10	<b>ESAOTE BIOMEDICA</b>								
32969	B800	roll	80	20	-	no	or/red	10	32980	P 80 base	roll	90	28	-	no	or/red	10
32951	BTL 08/MT, 08/MD, 08/MD3	roll	112	25	-	orange	orange	10	32981	P 80 base	pack	90	70	400	yes	or/red	25
32996	BTL 08/MT, 08/MD, 08/MD3	roll	58	25	-	no	or/red	20	32982	P 80 Power, P 80 SIX	pack	210	280	215	yes	or/red	1
32953	08 MT, MD	roll	112	23	-	no	or/red	10	32983	P8000 base	pack	80	70	300	yes	or/red	20
32952	BTL 08/LC, 08/LT	roll	210	25	-	no	or/red	5	32984	P8000 Power	pack	210	280	200	yes	or/red	1
32954	Atria 3100/6100, quest	pack	210	300	200	yes	or/red	1	32985	Archimed 4240, Archiwin	pack	210	300	200	yes	or/red	1
32955	Elan 100	pack	112	100	300	yes	or/red	25	32987	My Cardiopad base	pack	114	90	200	yes	or/red	20
32956	Delta 1 plus	pack	60	100	300	yes	or/red	25	33040	Esaoite My Cardiopad 2157036	pack	210	140	180	yes	or/red	5
32957	Delta 1 plus	roll	60	30	-	no	or/red	20	<b>FUKUDA DENSHI</b>								
32958	Delta 3 plus	pack	112	100	300	yes	or/red	25	32988	501A, FX121, FCP-500, OP18TE	roll	50	30	-	no	or/red	20
32960	Delta 60 plus, ECG 200L	roll	210	30	-	no	or/red	5	32989	FCP15, FCP11, OP110TE, OP119TE	roll	63	30	-	no	or/red	20
32961	Elan 2100, Delta 612	pack	210	280	200	yes	or/red	1	32994	FCP2101, 2155, FX 2111, CASE10, FCP 7101, FX 7102	roll	145	30	-	no	or/red	10
33053	ECG 100+, ECG 100S	pack	100	150	200	yes	or/red	10	32990	OP 111 TE, OP 115 TE, FX102	pack	50	75	400	yes	or/red	20
33052	ECG 100/L	roll	100	20	-	orange	orange	10	32991	OP 112 TE, FX2111, FX7101, FX7102, FCP2201 G	pack	63	75	400	yes	or/red	25
33054	AR2100 (VIEW/ADV) - 66010042	pack	210	140	200	yes	or/red	8	32992	FCP231, FCP4101, FX4010 FX7402, FX7542, FX8322, FDX4510	roll	210	30	-	no	or/red	5
33055	ECG 200+	pack	210	295	180	yes	or/red	1	32965	OP222TE, FX7202	pack	110	140	143	yes	or/red	20
32963	ECG-1101B, 1101C, 1101G	roll	50	25	-	no	or/red	20	<b>GIMA</b>								
33025	ECG-1103B, 1103G, 1103L	roll	63	30	-	no	or/red	20	32996	Cardiogima 1M, 3M	roll	58	25	-	no	or/red	20
33058	ECG-1103GW, 1103LW	roll	80	30	-	no	or/red	10	32953	Cardiogima 6M, colour	roll	112	23	-	no	or/red	10
33059	ECG-1106G, 1106L	roll	112	25	-	no	or/red	10	32950	Cardiogima 12, 12M	roll	215	25	-	no	or/red	5
32950	ECG-1112L, 112M, T/12	roll	215	25	-	no	or/red	5	<b>HEWLETT PACKARD</b>								
32963	CM100	roll	50	25	-	orange	orange	20	32999	Page Writer 100/200, M1709A, M2483A	pack	210	300	200	hole	or/red	1
32969	CM300, H3	roll	80	20	-	no	or/red	10	32988	40453 A, 43120	roll	50	30	-	no	or/red	20
32965	CM600	roll	110	20	-	no	or/red	20	<b>INNOMED</b>								
32992	CM 1200 A/B	roll	210	30	-	no	or/red	5	32996	HearthMirror IKO 3, HeartScreen 60 IKO, 60 G	roll	58	25	-	no	or/red	20
32966	CM 1200 A/B	pack	210	295	150	yes	or/red	1	32964	HeartScreen IKO 80G-L, 80C-L	roll	80	25	-	no	or/red	10
32963	CMS80 Cardiopoket	roll	50	25	-	orange	orange	20	32953	HeartScreen 112 C-D	roll	112	23	-	no	or/red	10
32968	90A, 100G	roll	50	20	-	no	or/red	20	33004	HeartScreen 210	pack	210	150	200	orange	orange	10
32969	300G, 300GA, 300GT	roll	80	20	-	no	or/red	10	33059	ECG-1106G - ECG-1106L	roll	112	25	-	orange	orange	10
32970	600G	roll	110	20	-	no	or/red	10	<b>KENZ - SUZUKEN</b>								
32967	1200, 1201, 1212G	roll	210	30	-	no	or/red	5	33002	Cardico 302/601	roll	112	27	-	no	or/red	10
32963	LTD400	roll	50	25	-	no	or/red	20	33003	Cardico 302/601	pack	112	90	200	yes	or/red	20
32996	LTD410, HeartScreen 60G	roll	58	25	-	no	or/red	20	33058	Cardico 306	roll	80	30	-	orange	orange	10
32964	LTD405S, LTD405W, LTD415 Heartscr. 80G	roll	80	25	-	no	or/red	10	32960	Kenz-Cardico 1210, 1211	roll	210	30	-	no	or/red	5
32951	LTD418, HeartScreen 112	roll	112	25	-	no	or/red	10	33025	ECG 305	roll	63	30	-	orange	orange	20
32965	LTD440, LTD445/S	pack	110	140	143	yes	or/red	20	32988	103, R50-30R, TEC 8300	roll	50	30	-	no	or/red	20
32960	LTD450, LTD455/S, LTD455/W	roll	210	30	-	no	or/red	5	<b>LUMED</b>								
									33058	EUROECG 3view	roll	80	30	-	orange	orange	10
									<b>MARQUETTE - HELIGE - GE</b>								
									33010	MAC PC	roll	108	23	-	no	or/red	10
									33011	Case 10-12-15, MAC 3000-5000	pack	215	280	300	hole	or/red	1
									32966	Cardiosmart, MAC 1200/1600	pack	210	295	150	yes	or/red	1

\*or/red: grid could be orange or red colour