

SPECIFICATION

Topper 8 3M non sterile



Universal swabs, non-sterile

4-ply, cross-folded. Made of 65% pure
cellulose viscose and 35% polyester.

Bleached without chlorine.

Satisfies DIN EN 1644-1:1997-04.

SPECIFICATION**Topper 8 3M non sterile**

Dimensions	OC	PZN	REF	Unit
5 x 5 cm	50 x 100 pcs.	03420464	J7 50405	100 pcs.
7,5 x 7,5 cm	56 x 100 pcs.	03420470	J7 50407	100 pcs.
10 x 10 cm	40 x 100 pcs.	03420487	J7 50410	100 pcs.
5 x 5 cm	42 x 200 pcs.	03420292	J7 53405N	200 pcs.
7,5 x 7,5 cm	28 x 200 pcs.	03420300	J7 53407	200 pcs.
10 x 10 cm	20 x 200 pcs.	03420317	J7 53410	200 pcs.
10 x 20 cm	10 x 200 pcs.	03420323	J7 53420	200 pcs.

