

## BATH/SHOWER CHAIRS



• **27730 SHOWER CHAIR with backrest**  
Aluminium frame, seat and backrest in plastic material. Height adjustable (35-45 cm), non slip end. Load: 100 kg.



• **42906 REINFORCED SHOWER CHAIR with backrest**  
Aluminium frame, seat and backrest in plastic material. Height adjustable (42-52 cm), non slip end. Load: 150 kg.



• **42907 SHOWER CHAIR with backrest and armrest**  
Aluminium frame, seat and backrest in plastic material. Height adjustable (41-54 cm). Load: 100 kg.



Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR.  
• **42909 SHOWER CHAIR with PU backrest and seat**  
Aluminium silver frame with PU backrest and seat. Height adjustable foot tube. Load: 136 kg. Weight: 4 kg.

## SHOWER COMMUNE CHAIR

• **42913 AQ-TICA SHOWER AND TOILET CHAIR - dark grey/white**  
Two function chair: body care and personal hygiene. It can be used as a toilet chair or as a shower chair. Delivered as a kit in a retail box can easily be mounted without any tools. Rust free structure made of aluminium and polypropylene. Supplied with hygienic pan with lid.



## ROTATING SHOWER CHAIR

• **42912 ROTATING SHOWER CHAIR**  
Height adjustable and 360° rotating seat chair. Swing-up armrests make patient transfer easier. The lever locks the seat in 4 different positions. Aluminium silver frame with comfortable PU padded seat and backrest. Foam arm pad. It can be folded.  
Seat dimension: 40.7x33.7x h 5 cm  
Max load: 136 kg. Weight: 5.3 kg.  
Packing size: 52.5x36x h 42.5 cm  
Folding dimension: 51x34.5xh 40 cm.  
Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR.




## SHOWER COMMUNE WHEELCHAIR

• **42915 SHOWER COMMUNE WHEELCHAIR**  
Foldable shower commode wheelchair with white painted aluminium frame and plastic seat. Flip-up armrest and lockable front castors. Size: 75x44x48 cm  
Packing size: 67x29x53 cm  
Weight: 6 kg.  
Load: 100 kg.



## BATH CHAIR

• **27755 SWIVEL BATH CHAIR**  
Revolving bath chair at 360° with 4 positions for easier access to bath. Polypropylene seat. Frame in painted steel, rotation system in stainless steel. Load: 100 kg. Weight: 5.5 kg.



## BATH SEAT

• **42910 BATH SEAT - adjustable**  
Aluminium frame and plastic seat (41x23 cm)  
Load: 100 kg.

