## **MEDICAL SCALES 2024**

**Professional Care** 













### Convenient circumference measuring tape with practical pull-out mechanism

### Features

#### **Technical data**

KERN

• Continuous pull-out measuring tape for precise measurement e.g. around the body

000

- Ideal to determine the WHR (Waist-to-Hip-Ratio), which can give information about the distribution of the body's fat deposits using an indirect determination of abdominal fat
- · With integrated feed mechanism
- The easiest handling:
  - Pull out the measuring tape
- Hang it in the recess provided for this
- Operate the feed mechanism
- Read measurement value with millimetre precision
- Note: The models are only delivered in a set of 5 units. That means the price refers to a delivery of 5 units. Cannot be delivered individually

• Measuring tape length 2,05 m (80 inch)

• Net weight approx. 0,25 kg (per piece)


STANDARD				
1 DAY				
1 DAY				

Model	Measuring range	Increments of	Overall dimensions	Net weight approx.	
			W×D×H		
KERN	cm (inch)	mm	mm	kg	
MSW 200S05	205 (80)	1	55×90×20	0,25	

# **MEDICAL SCALES 2024**

#### **KERN Pictograms**





one eye

With 3 W LED illumination Darkfield condenser/unit For a higher contrast due to Infinity corrected optical Automatic temperature

For measurements between

conformity assessment is specified in the pictogram

the pictogram Pallet shipment

The time required for internal shipping preparations is shown in days in the pictogram

\*The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by KERN & SOHN GmbH is under license

Other trademarks and trade names are those of their respective owner

