

KERN MGD 200K-1L

KERN

Elegant look combined with large, flat glass weighing plate to make it easier to step on and give maximum safety



Construction

Dimensions weighing platform (W×D×H)	380×330×27 mm
Dimensions completely mounted (W×D×H)	380×330×27 mm
Dimensions weighing surface (W×D)	380×330 mm

Functions

Auto-off interval(s) in battery mode/rechargeable battery mode	1 min
--	-------

Power Supply

Supplied power supply	Battery
Battery	3×1,5V AAA
Battery / accumulator type	Alkali(-Manganese)
Battery connection	Battery insert
Battery capacity	420 mAh
Battery voltage	1,5 V
Battery operating time	100 h

Environmental conditions

Humidity of environment [Min]	5 %
Humidity of environment [Max]	80 %
Ambient temperature [Min]	-10 °C
Ambient temperature [Max]	40 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	410×352×43 mm
Net weight	1,5 kg
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	2,8 kg
Shipping weight	2,8 kg

Services (optional)

Article number for DAkkS calibration	963-129
Article number for certificate of conformity (verification)	969-517

Category

Brand	KERN
Product category	Medical scale
Product group	Bathroom scale
Product family	MGD

Measuring System

Construction type of the scale	Single-range balance
Weighing system	Strain gauge
Weighing capacity [Max]	250 kg
Readability [d]	0,1 kg
Resolution	2.500
Linearity	± 1,35 kg
Reproducibility	1,35 kg
Default unit	kg
Units	kg lb
Adjustment options	Calibration not possible
Stabilization time	2 s
Warm-up time	10 min
Eccentric loading at 1/3 [Max]	1,35 kg
Maximum creep (15 minutes)	200 g
Maximum creep (30 minutes)	400 g
Tolerance (% of [Max])	0,5%

Approval

CE mark	✓
---------	---

Display

Display type	LCD
Display digit height - smallest digit	10 mm
Display digit height	50 mm

KERN MGD 200K-1L



Elegant look combined with large, flat glass weighing plate to make it easier to step on and give maximum safety

Pictograms

STANDARD



OPTION

