

KERN MSF 200

KERN

When every centimetre counts – mechanical height rod approved as a medical device for professional use in medical diagnostics



Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	1015×70×45 mm
Net weight	0,789 kg
Shipping method	Parcel service
Net weight approx.	0,80 kg
Gross weight approx.	1,2 kg
Shipping weight	1,006 kg

Pictograms

STANDARD



Category

Brand	KERN
Product category	Medical scale
Product group	Height rod
Product family	MSF

Measuring System

Default unit	mm
Units	inch/ft cm
Division	0,1 cm
Measuring range length [Min]	60 cm
Measuring range length [Max]	205 cm
Tolerance measuring rod	10 mm
Measuring range measuring rod (cm)	60-100cm @ 0,1cm ± 0,5cm 100-210cm @ 0,1cm ± 1cm

Approval

CE mark	✓
Approved as a medical device	✓

Construction

Dimensions device (W×D×H)	60×60×930 mm
Dimension housing (W×D×H)	53×40×1000 mm
Measurement alignment	vertical

Environmental conditions

Use temperature [Min]	10 °C
Use temperature [Max]	40 °C