

# KERN MSW 200

**KERN**

Convenient circumference measuring tape with practical pull-out mechanism



Net weight approx.	0,05 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Medical scale
Product group	Body size measurement
Product family	MSW

## Measuring System

Default unit	mm
Units	mm
Division	0,1 cm
Measuring range length [Min]	10 cm
Measuring range length [Max]	205 cm
Tolerance measuring rod	5 mm
Measuring range measuring rod (cm)	0-205cm @ 0,1cm ± 0,5cm

## Construction

Dimensions device (W×D×H)	55×90×20 mm
Dimension housing (W×D×H)	55×90×20 mm
Measurement alignment	horizontal

## Environmental conditions

Use temperature [Min]	10 °C
Use temperature [Max]	40 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	70×100×25 mm
Net weight	0,045 kg
Shipping method	Parcel service